














Sound Delivery Compatibility Guide

	 Audéo MINI SDS 1.5	 Audéo MINI SDS 2.0	 Audéo SMART SDS 1.5	 Audéo SMART SDS 2.0	 Audéo YES SDS 1.5	 Audéo YES SDS 2.0	 CORE micro*
 Original SlimTip <i>Soft or Hard</i>		Standard Receiver		Standard Receiver	Standard Receiver		■
 SlimTip Extra Retention <i>Soft or Hard</i>		Standard Receiver		Standard Receiver	Standard Receiver		■
 cShell		Standard Receiver Power Receiver		Standard Receiver Power Receiver	Standard Receiver Power Receiver		
 Closed Dome		Standard Receiver Power Receiver		Standard Receiver Power Receiver	Standard Receiver Power Receiver		■
 Open Dome		Standard Receiver		Standard Receiver	Standard Receiver		■
 Power Dome		Standard Receiver Power Receiver		Standard Receiver Power Receiver	Standard Receiver Power Receiver		■

■ compatible

*Exélia Art, Versáta, Certéna, Nios, Milo Plus



xReceiver Standard



xReceiver Power



SDS 1.5 Connector



SDS 2.0 Connector



Audéo Family Quick Reference Guide

Feature/Benefit	Premium			Advanced			Basic		
	Audéo MINI IX	Audéo SMART IX	Audéo YES IX	Audéo MINI V	Audéo SMART V	Audéo YES V	Audéo MINI III	Audéo SMART III	Audéo YES III
SoundRecover improves audibility by compressing high frequency sounds	■	■	■	■	■	■	■	■	■
SoundFlow seamlessly and automatically adapts in changing environments	SoundFlow Premium 4 base programs			SoundFlow Advanced 3 base programs			SoundFlow Standard 2 base programs		
WhistleBlock provides feedback free hearing while providing more dynamic stable gain across a wider frequency range	■	■	■	■	■	■	■	■	■
VoiceZoom increases speech understanding in noise with multiple channel directionality	33 channels			12 channels			Digital AudioZoom		
NoiseBlock provides premium comfort by automatically reducing background noise	20 channels			16 channels			6 channels		
Real Ear Sound restores front/back localization	■	■	■	■	■	■			
WindBlock enhances outdoor enjoyment by suppressing wind noise without effecting sound quality	■	■	■	■	■	■			
EchoBlock improves understanding in reverberant environments by reducing echo	■	■	■						
SoundRelax eliminates annoying sounds (clanging dishes)	■	■	■						
Acoustically Optimized Vent (with SlimTip or cShell) offers maximum comfort with minimum occlusion	■	■	■	■	■	■	■	■	■
QuickSync provides convenient one touch binaural control of volume and/or programs		■			■			■	
ZoomControl provides wearers the ability to choose their listening focus, ideal in a car		■	■ via myPilot						
DuoPhone automatically enhances phone signals in both ears with an optimal signal		■							