

Phonak CROS


Discover the unique solution for single-sided hearing from Phonak





Discover all around hearing delight

Phonak CROS is the smallest and most stylish wireless CROS/BiCROS solution designed to meet your individual needs and preferences. Phonak CROS consists of only two parts – a transmitter microphone for the ear that cannot benefit from a hearing instrument, and a receiver hearing instrument on the better hearing ear. This means sound from your poorer hearing side is wirelessly sent to your better ear.





CROS or BiCROS

When your better ear hears normally



CROS

A small transmitter, placed behind or inside the ear that cannot benefit from a hearing aid, picks up sound and transfers it wirelessly to a receiver on your normal hearing ear.

This will allow you to hear effortlessly even when sounds come from the side where you have no hearing.

When you have a hearing loss in your better ear



BiCROS

The small transmitter placed behind or inside the ear that cannot benefit from a hearing aid picks up sound and transfers it wirelessly to your ear fitted with a

hearing aid. At the same time this hearing aid amplifies sounds coming directly to this ear. Your hearing and understanding experience will improve and your communication capabilities will be better than ever.



Phonak CROS transmitter styles



Phonak CROS H20

Together with a WaterResistant Phonak hearing aid, this robust little system gives you the freedom to enjoy your favorite activities with great confidence. Be it around water, on the beach, working out, in the rain or if you simply forget to take your system off to take a shower, there is no need to worry.



Phonak CROS

If you appreciate a small, stylish design with all the benefits that Phonak CROS can provide, this is an attractive solution for you.



Phonak CROS 13 and CROS 312

For a custom solution, a Phonak CROS in-the-ear option is for you. These custom-made devices come with a choice of battery sizes so that you can choose to have the benefit of a smaller device or a longer lasting battery.

A truly delightful hearing experience

"My daughter Halie recently got the Phonak CROS and it has truly been life changing! For the first time she feels as if she is a 'normal' hearing child. She has always had a severe hearing loss and has been wearing hearing aids for some years now, but the left one never really did anything for her. The day she received them was amazing. Her doctor placed them in her ears and I then whispered in her left ear 'I love you'. She turned her head and said 'I love you too mom'.

This aid has given her her confidence back. Her dad will be returning from deployment in the next few weeks and will be able to watch our baby 'hear' for the first time. We are looking forward to ordering the FM system in the near future which will only make this dream even better."

Angie Kitchens – mother of Halie, Phonak CROS wearer



Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

Interact freely. Communicate with confidence.
Live without limit. Life is on.

www.phonak-us.com

