ReSound Alera®

More models to fit more patients.



NEW ReSound Alera Traditional Custom Products

ReSound Alera Custom Remote Microphone Products

ReSound Alera Power BTE

NEW ReSound Alera Mini BTE

NEW ReSound Alera RIE

ReSound Alera Mini RIE

ReSound Alera Micro RIE





ReSound Alera®

Introducing...

RIE, Mini BTE and Traditional Custom Products New models to fit more patients

ReSound Alera offers a full complement of hearing aid models, each offering unique benefits. Alera fits 98% of all hearing losses and provides the kind of accurate, personalized fittings that are fundamental to good hearing and excellent sound. Alera also addresses a wide range of personal preferences and gives every client a comfortable, easy-to-handle solution.

NEW! WindGuard™ Incredibly clean signal

WindGuard™ is a unique wind-noise algorithm that takes advantage of inputs from dual microphones to provide patients with effective wind noise reduction. WindGuard™ seamlessly transitions from one environment to the next to maximize comfort and prevent sudden volume changes as the patient moves in and out of windy conditions. All of this is done without compromising speech intelligibility.

NEW! In-situ Audiometry More convenient than ever

In-situ audiometry is a convenient addition to ReSound Alera's already robust capabilities. With ReSound Alera, you can obtain an audiogram using the convenient in-situ audiometry using the patient's own Alera BTE or RIE. In-situ audiograms are completely interchangeable with conventional ones.

Alera

ITE 50

ReSound Alera

ReSound

Alera

RM 10

Custom Remote Microphone Products

ReSound

Alera

RM 30

ReSound

Alera

RM 50

ReSound Alera lineup

Alera

MC 20

Alera

CIC 10

ReSound Alera Traditional Custom Products NEW NEW NEW ReSound ReSound ReSound ReSound

Alera

ITF 40

Alera


ITC 30


Facts

- iSolate[™] nanotech protection
- Fits mild to profound hearing losses
- Receiver-in-ear, Mini BTE and Power BTE options
- Unique custom products with natural wind noise protection, directivity and localization
- 2.4GHz wireless connectivity*

* Selected models

Want to know more? www.nostreamers.com







ReSound Alera

Transitional BTE

ReSound Alera

Standard and Power 77/87







ReSound Alera Mini BTF 67

ReSound Alera

NEW

Mini BTE

ReSound Alera 62

ReSound Alera

RIE

ReSound Alera 61

Alera 60