

ReSound LiNX™

Say hello to the world's smartest hearing aid



ReSound



red|discover hearing



Imagine a new kind of hearing aid that lets you hear more of everything.

With ReSound LiNX™ everything you hear is vivid and natural. Imagine speech that's clear, strong and easy to understand no matter where you are, what you're doing or who you're listening to.

Imagine a video call with a loved one on your iPhone®, watching a favorite video on your iPad® or listening to the music that moves you on your iPod touch® all through a wireless stereo headset that just happens to be your hearing aids.


Imagine a hearing aid personalized to your every need and so small no one ever needs to know you're wearing it.

Imagine all that and you've only just begun to imagine ReSound LiNX.

ReSound LiNX is the world's first hearing aid made for your iPhone, iPad and iPod touch.

Visit resoundlinx.com to experience the world's smartest hearing aid.





"brilliant"

It sounds so natural

In a recent independent study comparing premium products from six manufacturers, a panel of experienced hearing-aid wearers gave ReSound's sound quality top rating.*

We think you will too. We also think you'll be able to relax and enjoy yourself in situations where it used to be tough – and tiring – to follow the conversation.

*SenseLab hearing aid benchmark test, 2013.

Read more at resoundlinx.com



"no kidding"

A direct connection to everything you love

ReSound LiNX is the first hearing aid made to connect directly to your Apple iPhone, iPad and iPod touch.* So you can explore the endless possibilities of streaming telephone conversations, music and other high-quality sound directly to your hearing aids.



**ReSound LiNX is compatible with iPhone 5s, iPhone 5c, iPhone 5, iPad Air, iPad (4th generation), iPad mini with Retina display, iPad mini, and iPod touch (5th generation) using iOS 7.X or later.*

Read more at resoundlinx.com





"incredible"

Hear more than anyone

ReSound Unite™ wireless accessories let you hear things that even people with normal hearing sometimes have trouble catching. Like a quiet voice from across a crowded room or the radio in the kitchen. You can even chat hands-free on your smart phone with anyone. Anywhere. Anytime.



Read more at gnresound.com/unite

It will look great on you



An all-new design makes ReSound LiNX nearly invisible - unless you want to show it off. It's small, slim and hugs the contours of your ear – for style, comfort and stability.



ReSound LiNX is a mini receiver-in-the-ear (RIE) hearing aid. It's worn behind the ear with a thin tube that wraps discreetly along the front of the ear and ends with a tiny bud inside your ear.



ReSound LiNX comes in a choice of 10 attractive colors (not shown in actual size)

** The ReSound Smart app is compatible with iPhone 5s, iPhone 5c, iPhone 5, iPad Air, iPad (4th generation), iPad mini with Retina display, iPad mini, and iPod touch (5th generation), using iOS 7.X or later. If you have an iPhone 4S or earlier, or an Android smartphone, you can use the ReSound Control app together with the ReSound Unite Phone Clip+. See resound.com/control*

Personalized just for you

The ReSound Smart™ app

Your hearing care professional customizes your ReSound LiNX hearing aids to precisely match your hearing loss. And with the new ReSound Smart app, you can personalize and control your hearing aids from your phone.*



Read more at gnresound.com/smartapp

It won't cramp your style

You might expect a hearing aid this good looking to be delicate, but ReSound LiNX is coated with iSolate™ nanotech, an ultrathin polymer shield that seals off every component from moisture, oil and dust. So it's rugged and durable.



ReSound



ReSound® helps people rediscover the joy of hearing.
We develop hearing solutions that imitate the natural ear
enabling them to live rich, active and fulfilling lives.

resoundlinx.com

***ReSound LiNX products are only available through your
local, independent hearing healthcare professional.***

ReSound North America

8001 Bloomington Freeway
Bloomington, MN 55420
1-888-735-4327
customerexperience@gnresound.com
gnresound.com

ReSound Canada

303 Supertest Road
Toronto, Ontario M3J 2M4
1-888-737-6863
gnresound.ca

ReSound Government Services

8001 Bloomington Freeway
Bloomington, MN 55420
1-800-392-9932
VAinfo@gnresound.com

MK603991RevA 01.2014

©2014 The GN ReSound Group, all rights reserved. Apple, iPad, iPhone and iPod touch are trademarks of Apple Inc, registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android and Google Play are registered trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. Bluetooth is a trademark of Bluetooth SIG, Inc.